

**IT - NOTA INFORMATIVA GUANTI MONOUSO**  
**EN - INSTRUCTION FOR USE FOR SINGLE USE GLOVES**  
**ES - FOLLETO INFORMATIVO GUANTES DESECHABLES**  
**FR - MODE D'EMPLOY GANTS À USAGE UNIQUE**  
**PT - INSTRUÇÕES DE UTILIZAÇÃO LUVAS DESCARTÁVEIS**  
**RO - INSTRUCȚIUNI de UTILIZARE MĂNUȘI de UNICĂ FOLOSINȚĂ**  
**DE – GEBRAUCHSANLEITUNG EINMALHANDSCHUHE**

1	 Reg. UE 2016/425	2	 EN 16350:2014	3	 1935/2004 – 2023/2006
4	 EN ISO 21420:2020	5	 EN 421:2010	6	 ISO 18889:2019
7	 A B C D E P EN 388:2016+A1:2018	8	 Type A, B or C EN ISO 374-1:2016	9	 ( VIRUS ) EN ISO 374-5:2016
10		11		12	

Art.	Materiale / Material / Matériau /	Finitura / Finish / Acabado / Finition	Interno / Interior / Interioer
GGU501	Lattice Naturale / Natural Latex / Látex natural / Latex Naturel / Naturlatex	Testurizzata / Textured / Texturé Texturizado / Textură / Texturiert	Clorinato / Chlorinated Clorado / Chloré / Clorurat Gechlort
GGU502	Nitrile / Nitrile / Nitriilo / Nitrile / Nitril	Diamante / Diamond / Diamant	
GGU508			

13	Modulo B / Module B Módulo B / Modulul B / Modul B	Numero / Number Número / Nombre Număr / Nummer	Indirizzo / Address / Dirección / Adresse Endereço / Adresa / Adresse
	SATRA Technology Europe Limited	2777	Bracetown Business Park. Clonee. D15Y2NP – Ireland

14	Modulo C2 / Module C2 Módulo C2 / Modulul C2 / Modul C2	Numero / Number Número / Nombre Număr / Nummer	Indirizzo / Address / Dirección / Adresse Endereço / Adresa / Adresse
	SATRA Technology Europe Limited	2777	Bracetown Business Park. Clonee. D15Y2NP – Ireland

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**EN** Gloves are designed as a single use two way barrier protection against cross contamination and also protect the hands mainly against chemical splash risks and comply with the standards shown by the relevant pictograms. Multipurpose: suitable for industrial activities, for laboratory use, quality control, gardening, DIY and hobby working. – **USE:** this Instructions for Use must be used in combination with the specific information that appears on the dispenser. These gloves are designed to protect the hands against the risks as shown by the pictograms depicted, as defined in the relevant EN or ISO standards. Ensure that the products are used only for their designated purpose. – **MARKINGS in ACCORDANCE WITH LAW:** (1) **CE:** the glove meets the essential health and safety requirements of EU Regulation 2016/425. Subjected to and complies with an assessment in accordance with the harmonized technical standards listed below. – **PPE Cat III:** Personal Protective Equipment of complex design that protects the user against fatal and serious risks that are. EU Type Examination Certificate (module B) issued by (13). Conformity to type based on quality assurance of the production process (Module D) issued by (14). Marks (1, 2) are followed by a four digit code which refers to the identification number of the Notified/Approved Body that is in charge of the PPE category III conformity assessment. The EU or UK Declaration of Conformity can be found on [rossini1969.it](http://rossini1969.it) – **EXPLANATION of MARKINGS and PICTOGRAMS** that may appear on gloves/packaging: product complies with Technical Standard (4) **EN ISO 21420:2020** [general requirements protective gloves]: read the instructions for Use prior to using the products, or contact ROSSINI for further information. If level "X" is mentioned under any of the pictograms, this means this test is not applicable and glove is not designed and therefore not be used for the specific risk. Level "0" indicates that the glove's performance is below the minimum level for a given risk. – (7) **EN 388:2016+A1:2018** [protection against mechanical hazards]: (A) Abrasion resistance (performance levels 0 to 4) / (B) Blade cut resistance (performance levels 0 to 5) / (C) Tear resistance (performance levels 0 to 4) / (D) Puncture resistance (performance levels 0 to 4) / (E) TDM EN ISO 13997 cut resistance (performance levels A to F) / (P) Impact protection (optional): impact protection in the knuckle area (not applicable to the finger area, which cannot be tested). If no P is declared, no impact protection is provided. **WARNING!** the performance (A to E) claimed for gloves are based on tests performed on the palm area of the gloves only. – (9) **EN ISO 374-5: 2016** [protection against bacteria, fungi and virus]. – (8) **EN ISO 374-1:2016 TIPO A, B or C** [protection against chemical hazards]: **Type A:** chemical breakthrough time > 30 minutes against at least six chemicals as per list defined in EN ISO 374-1:2016. **Type B:** chemical breakthrough time > 30 minutes against at least three chemicals as per list defined in EN ISO 374-1:2016. **Type C:** chemical breakthrough time > 10 minutes against at least one test chemical as per list defined in EN ISO 374-1:2016 (no code underneath the pictogram). A = methanol, B = acetone, C = acetonitrile, D = dichloromethane, E = carbon disulfide, F = toluene, G = diethylamine, H = tetrahydrofuran, I = ethyl acetate, J = n-heptane, K = sodium hydroxide, 40%, L = sulphuric acid, 96 %, M = nitric acid, 65%, N = acetic acid, 99%, O = ammonia, 25%, P = hydrogen peroxide, 30%, S = hydrofluoric acid, 40%, T = formaldehyde, 37%. – (6) **ISO 18889: 2019** [protection against pesticides]: suitable when the potential risk is relatively low. These gloves are not suitable for use with concentrated pesticide formulations and/or for scenarios where mechanical risks exist. The pesticide shall not have the possibility to penetrate between the garment sleeve and the glove. If the overlap is less than approximately 50 mm between the glove and the sleeve, a glove with a longer length should be used. For pesticide resistance data, the duration of the test is not based on actual use time since the permeation test is an accelerated test in which the surface of the specimen is in constant contact with the test chemical. Although the duration of the exposure may be for a longer period during field application with a dilute formulation, the entire surface is not in constant contact with the test chemical. **WARNING!** Chemical resistance data has been assessed under laboratory conditions and relates only to the chemical tested. It can be different if used in a mixture. For gloves equal or longer than 400 mm, the chemical resistance data is based from samples taken 80 mm from the end of the cuff. The data may not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. Check that the gloves are suitable for the intended use. Conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Chemical permeation data and degradation data are available upon request. – (5) **EN 421:2010** [protection against radioactive contamination]. – (3) **FOOD:** suitable for contact with foodstuffs. Products carrying this pictogram are in conformity with the European Regulations 1935/2004 and 2023/2006 as well as with all applicable National Regulations for Food-contact materials. – (12) **Non-Sterile**. – (11) **Made with / Without natural rubber latex**. – (10) **Single Use Only**. – (10) **Eprly and Manufacture date**. – **PRECAUTIONS FOR USE:** before usage, inspect the gloves for any defects or imperfections such as holes, pinholes and tears. If the gloves are ripped or punctured during use, dispose of them immediately. If in doubt, do not use the gloves, get a new pair. It is essential to keep all chemicals from contact with the skin, even if they are thought to be harmless. Ensure the chemicals cannot enter via the cuff. Remove the glove immediately if contaminated by a concentrated spill of pesticides. Contaminated gloves should be cleaned or washed or wiped dry before removal. Avoid touching contaminated surfaces with bare hands. – **Doffing:** slide each glove onto your hand, ensuring a comfortable fit without any wrinkles that could impede movement or grip. Adjust the gloves as needed to ensure they cover your hands fully and securely. **Doffing:** gloves must be cleaned, washed or wiped dry before removal. Pinch and grip with one hand the gloved fingers of the opposite hand. Slowly pull half glove off, being careful not to touch the interior surface or the hand with the dirty glove. With half glove off repeat the previous step for he opposite hand. Hold the removed gloves from the cuff. – Avoid touching contaminated surfaces with bare hands. Gloves should not come in contact with a naked flame nor used for protection against heat. Gloves shall not be used for protection against ionizing radiation nor for use in containment enclosures. Gloves suitable for contact with foodstuffs may show some migration against specific foodstuffs. If gloves are being used in explosive environments (ATEX), please ensure they meet the (2) **EN 16350:2014** [electrostatic properties] requirements. Persons wearing these gloves should be properly earthed, e.g. by wearing adequate footwear & clothing. **WARNING!** the gloves shall not be unpacked, opened, adjusted or removed whilst in flammable or explosive atmospheres. The electrostatic properties of the gloves might be adversely affected by ageing, wear, contamination and damage and might not be sufficient for oxygen enriched flammable atmospheres where additional assessments are necessary. – **COMPOSITION/HAZARDOUS INGREDIENTS:** some gloves may contain components known to be potentially allergenic in individuals sensitized to such components. These persons may develop irritation and/or allergic contact reactions including anaphylactic responses. If allergic reaction should occur, obtain medical advice. **WARNING!** If the glove contains natural latex the symbol (11) is on the packaging, these gloves MAY CAUSE ALLERGIC REACTIONS in sensitized individuals. – **CARE INSTRUCTIONS. Storage:** do not expose to direct sunlight; store in a cool, dry place in the original packaging. Keep away from heat and ozone sources. If properly stored, the gloves will not lose their performances or undergo any significant change in their characteristics for up to 3 (three) years from the date of manufacture. **Cleaning:** DO NOT wash; DO NOT re-use. **Disposal:** used gloves contaminated with infectious or other hazardous materials must be disposed and not reused. Gloves must be disposed once they show visible signs of degradation during use, such as discoloration, tearing, holes and weakening. Dispose in accordance with local regulations. Dispose in a landfill site or incinerate under controlled conditions. – **NOTE:** the information contained in this document together with the results obtained from the physical examination performed by the Notified Body should assist in glove selection. However, they do not reflect the actual protection of gloves in the workplace due to other factors such as temperature, abrasion, degradation, etc., in accordance with EU regulation

2016/425. These products have been manufactured under a registered quality system and comply with the requirements established by the ISO 9001:2015 standard. None of the materials or processes used to manufacture these gloves are known to be harmful to the user.

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